

**T4 DNA Polymerase (5 U/μL)****Thermo Scientific**

T4 DNA polimerasi è template-dipendente e catalizza la sintesi 5'-3' da singoli filamenti di DNA innescati.



Code	Description	Packaging
FL5EP0061	T4 DNA Polymerase (5 U/μL)	100 U (*)
FL5EP0062	T4 DNA Polymerase (5 U/μL)	500 U (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

