



**RevertAid H Minus Reverse Transcriptase (200 U/ $\mu$ L)**

**Thermo Scientific**

M-MuLV RT ricombinante senza attività di RNasi H alta resa nella sintesi di cDNA.



| Code      | Description                                        | Packaging   |
|-----------|----------------------------------------------------|-------------|
| FL5EP0451 | RevertAid H Minus Reverse Transcriptase (200 U/μL) | 10000 U (*) |
| FL5EP0452 | RevertAid H Minus Reverse Transcriptase (200 U/μL) | 50000 U (*) |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

