



T4 RNA ligase

Thermo Scientific

T4 RNA Ligasi - catalizza la formazione ATP-dipendente di legami fosfodiesteri tra 5'-P e 3'-OH terminale di oligonucleotidi, RNA e DNA a singolo filamento.



Code	Description	Packaging
FL5EL0021	T4 RNA ligase	1000 U (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

