


**TRIPZ - Dharmacon# part of horizon#****Dharmacon# part of GE
Healthcare**

TRIPZ - vettori per il silenziamento stabile attraverso short harpin RNA (shRNA) inducibile - tecnologia microRNA-adapted - linea completa di genomi per uomo - la presenza del gene RFP permette di controllare l'efficienza dell'esperimento - in stock di glicerolo - utilizza plasmidi di origine lentivirale, l'utilizzo di particelle lentivirali necessita di un laboratorio con livello di biosicurezza S2.



Code	Description	Packaging
FE6RHS4743	TRIPZ Inducible Lentiviral Non-silencing shRNA Control (in glicerolo)	1 each (*)
FE6RHS4750	TRIPZ Inducible Lentiviral Empty Vector shRNA Control (in glicerolo)	1 each (*)
FE6RHS4741	TRIPZ Inducible shRNA Transfection Starter Kit (in glicerolo)	1 each (*)
FE6RHS5087	TRIPZ shRNA Packaging Starter Kit (in glicerolo)	1 each (*)
FE6RHS4696	TRIPZ Lentiviral shRNA Individual - uomo (in glicerolo)	1 each (*)
FE6RHS4744	TRIPZ uomo GAPDH Inducible Lentiviral shRNA Positive Control (in glicerolo)	1 each (*)
FE6RHS6371	TRIPZ Human Genome Inducible Lentiviral shRNA Library V.2	1 each



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

