


**TRC (pLOK1) - Dharmacon# part of horizon#****Dharmacon# part of GE
Healthcare**

TRC (pLOK1) - vettori per il silenziamento stabile attraverso short harpin RNA (shRNA) - linea completa per genomi di uomo e topo - in stock di glicerolo - utilizza plasmidi di origine lentivirale, l'utilizzo di particelle lentivirali necessita di un laboratorio con livello di biosicurezza S2.



Code	Description	Packaging
FE6RHS4875	Human Apoptosis TRC lentiviral shRNA gene family	1 each (*)
FE6RMM4954	TRC Apoptosis shRNA Gene Family Library (in glicerolo) - topo	1 each (*)
FE6RHS3979	TRC Lentiviral shRNA Individual (in glicerolo) - uomo	1 each (*)
FE6RHS4012	TRC Lentiviral Human Genome shRNA Library	glycerol stock
FE6RHS4332	TRC Lentiviral Druggable Gene Family Library (in glicerolo) - uomo	1 each (*)
FE6RHS4459	TRC Lentiviral eGFP shRNA positive control (in glicerolo)	1 each (*)
FE6RHS4533	TRC Lentiviral shRNA Target Gene Set (in glicerolo) - uomo	1 each (*)
FE6RHS4872	TRC Angiogenesis shRNA Gene Family Library (in glicerolo) - uomo	1 each (*)
FE6RHS4878	TRC Cell Cycle shRNA Gene Family Library (in glicerolo) - uomo	1 each (*)
FE6RHS4881	TRC GPCR shRNA Gene Family Library (in glicerolo) - uomo	1 each (*)
FE6RHS4884	TRC Kinase shRNA Gene Family Library (in glicerolo) - uomo	1 each (*)
FE6RHS4887	TRC Phosphatase shRNA Gene Family Library (in glicerolo) - uomo	1 each (*)
FE6RHS4890	TRC fattori di trascrizione shRNA Gene Family Library (in glicerolo) - uomo	1 each (*)
FE6RHS4893	TRC Ubiquitin shRNA Gene Family Library (in glicerolo) - uomo	1 each (*)
FE6RHS4896	TRC Ubiquitin Ligases shRNA Gene Family Library - (in glicerolo) - uomo	1 each (*)
FE6RMM3981	TRC Lentiviral shRNA Individual (in glicerolo) - topo	1 each (*)
FE6RMM4013	TRC Lentiviral Genome shRNA Library (in glicerolo) - topo	1 each (*)
FE6RMM4534	TRC Lentiviral shRNA Target Gene Set (in glicerolo) - topo	1 each (*)
FE6RMM4953	TRC Angiogenesis shRNA Gene Family Library (in glicerolo) - topo	1 each (*)
FE6RMM4955	TRC Cell Cycle shRNA Gene Family Library (in glicerolo) - topo	1 each (*)



Code	Description	Packaging
FE6RMM4956	TRC GPCR shRNA Gene Family Library (in glicerolo) - topo	1 each (*)
FE6RMM4957	TRC Kinase shRNA Gene Family Library (in glicerolo) - topo	1 each (*)
FE6RMM4958	TRC Phosphatase shRNA Gene Family Library (in glicerolo) - topo	1 each (*)
FE6RMM4959	TRC fattori di trascrizione shRNA Gene Family Library (in glicerolo) - topo	1 each (*)
FE6RMM4960	TRC Ubiquitin shRNA Gene Family Library (in glicerolo) - topo	1 each (*)
FE6RMM4961	TRC Ubiquitin Ligases shRNA Gene Family Library (in glicerolo) - topo	1 each (*)
FE6RMM5146	TRC Cancer Genes shRNA Library (in glicerolo) - topo	1 each (*)



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

