


**NucleoMag 96 Trace****Macherey Nagel**

NucleoMag 96 Trace - isolamento di DNA genomico da campioni forensi a scopo stoccaggio in database (prevalentemente da tampone buccale) - tecnologia: biglie magnetiche - formato: piastra a 96 pozzetti - rischio di contaminazione crociata minimizzato - piccoli volumi di eluizione: 50-200 μL - facilmente adattabile a sistemi automatici - resa dipendente dalla quantità e qualità del campione - concentrazione tipica: 10-30 $\text{ng}/\mu\text{L}$ - A260/A280 1,6-1,9 - tempo di preparazione: 120 min/preps - capacità legante: 0,4 $\mu\text{g}/\mu\text{L}$ biglie - comprende: piastra, biglie, tamponi, proteinasi K.



Code	Description	Packaging
FC144600C	Kit NucleoMag 96 Trace - con piastra a 96 pozzetti	4X96 PREP
FC144600G	Kit NucleoMag 96 Trace - con piastra a 96 pozzetti	24X96
FC1744600	Kit NucleoMag 96 Trace - con piastra a 96 pozzetti	96 prep

Accessories

Code	Description	Packaging
FC140481G	Square Well Block	24 pezzi
FC1740477	MN TUBE STRIP WITH CAP STRIP	4 set
FC1740481	SQUARE WELL BLOCK	4 pezzi



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

