

**NucleoSpin Tissue XS****Macherey Nagel**

NucleoSpin Tissue XS - purificazione di DNA genomico, batterico e virale da campioni piccolissimi - ideale per: analisi forense, materiale da biopsie, cellule da FACS - tecnologia: membrana di silice - formato: mini colonna, XS design: campioni: 0.025-10 mg di tessuti umani o animali, 1-30 μ L di sangue fresco o congelato, 10 10exp4 cellule umani o animali, 15-30 mm² di Guthrie card - dimensioni frammento: 20 bp- circa 50 kbo - volume di eluizione: 5-30 μ L - tempo di preparazione: 20 min/prep (esclusa lisi) - capacità legante: 50 μ g - comprende: colonne NucleoSpin® Tissue XS, provette per campionamento (2 mL), tamponi, proteinasi K.



Code	Description	Packaging
FC140901E	Kit NucleoSpin Tissue XS	10 PREP
FC140901L	Kit NucleoSpin Tissue XS	50 PREP
FC140901N	Kit NucleoSpin Tissue XS	250 PREP



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

