

**NucleoBond CB/AXG****Macherey Nagel**

NucleoBond CB/AXG - per la purificazione di DNA genomico da sangue intero - tecnologia: cromatografia a scambio anionico - formato: colonna - permette di ottenere DNA genomico ultra-puro a partire da differenti tipologie di campione (sangue intero, siero, plasma, colture cellulari, buffy coat) - dimensione dei frammenti: da 500 bp a 300 kbp - A260/A280 1,8/1,95 - tempo di preparazione: 4-5 h - dimensione dei campioni: fino a 20 mL - versione kit CB contiene: colonne NucleoBond AXG, tamponi, proteinasi K.



Code	Description	Packaging
FC1740508	Kit NucleoBond CB 100	20 prep.
FC1740544	NucleoBond AXG 20	20 pz
FC1740545	NucleoBond AXG 100	20 pz
FC1740546	NucleoBond AXG 500	20 pz

Accessories

Code	Description	Packaging
FC1740603	NUCLEOBOND BUFFER SET III	1 set
FC1740604	NUCLEOBOND BUFFER SET IV	1 set



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*