

**NucleoBond CB/AXG****Macherey Nagel**

NucleoBond CB/AXG - per la purificazione di DNA genomico da sangue intero - tecnologia: cromatografia a scambio anionico - formato: colonna - permette di ottenere DNA genomico ultra-puro a partire da differenti tipologie di campione (sangue intero, siero, plasma, colture cellulari, buffy coat) - dimensione dei frammenti: da 500 bp a 300 kbp - A260/A280 1,8/1,95 - tempo di preparazione: 4-5 h - dimensione dei campioni: fino a 20 mL - versione kit CB contiene: colonne NucleoBond AXG, tamponi, proteinasi K.



| Code      | Description           | Packaging |
|-----------|-----------------------|-----------|
| FC1740508 | Kit NucleoBond CB 100 | 20 prep.  |
| FC1740544 | NucleoBond AXG 20     | 20 pz     |
| FC1740545 | NucleoBond AXG 100    | 20 pz     |
| FC1740546 | NucleoBond AXG 500    | 20 pz     |

## Accessories

| Code      | Description               | Packaging |
|-----------|---------------------------|-----------|
| FC1740603 | NUCLEOBOND BUFFER SET III | 1 set     |
| FC1740604 | NUCLEOBOND BUFFER SET IV  | 1 set     |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

