

**NucleoSpin Plasmid EasyPure****Macherey Nagel**

NucleoSpin Plasmid EasyPure - per la preparazione ultra-veloce, su bassa scala, di DNA plasmidico ultra-puro: fino a 35 µg - LyseControl per la verifica del completamento della lisi alcalina - Liquid RNase A: facile da maneggiare, senza necessità di dissoluzione - la combinazione degli step di lavaggio e essiccazione permette una procedura più veloce: 6 mini preps in 14 min - tecnologia: membrana di silice - formato: mini colonna - campione: 2-10 mL di coltura E. coli - dimensione del vettore: 15 kbp - resa tipica: 15-30 µg - A260/A280 1,8-1,85 - volume di eluizione 50 µL - comprende: colonne NucleoSpin Plasmid® EasyPure, provette di campionamento (2 mL), tanponi, RNase A liquida.



Code	Description	Packaging
FC140727L	Kit NucleoSpin Plasmid EasyPure	50 reaz.
FC140727N	Kit NucleoSpin Plasmid EasyPure	250 reaz.

Accessories

Code	Description	Packaging
FC1740953	NUCLEOSPIN BUFFER SET	1 set



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

