

**NucleoBond PC/AX****Macherey Nagel**

NucleoBond PC/AX - per la purificazione di DNA plasmidico per transfezione - tecnologia: cromatografia a scambio anionico - formato: colonne di diverse dimensioni - versione in kit NucleoBond PC contiene: colonne AX, tamponi, RNasi A.

In particolare:

NucleoBond® PC 20: colonne NucleoBond® AX 20, tamponi, RNase A.

NucleoBond® PC 100: colonne NucleoBond® AX 100, tamponi, RNase A, NucleoBond® Folded Filters.

NucleoBond® PC 500: colonne NucleoBond® AX 500, tamponi, RNase A, NucleoBond® Folded Filters XL.

NucleoBond® PC 2000: colonne NucleoBond® AX 2000, tamponi, RNase A, NucleoBond® Folded Filters XL.

NucleoBond® PC 10000: colonne NucleoBond® AX 10000, tamponi, RNase A, NucleoBond® Folded Filters (Type 1, Type 2)



Code	Description	Packaging
FC140521M	Kit NucleoBond AX 100	100 pz
FC140531L	Kit NucleoBond AX 500	50 pz
FC140571F	Kit NucleoBond PC20	20 prep
FC140571M	Kit NucleoBond PC20	100 prep
FC140573F	Kit NucleoBond PC100	20 prep
FC140573M	Kit NucleoBond PC100	100 prep
FC140574E	Kit NucleoBond PC500	10 prep.
FC140574H	Kit NucleoBond PC500	25 prep
FC140574L	Kit NucleoBond PC500	50 prep
FC140574M	Kit NucleoBond PC500	100 prep
FC1740511	NucleoBond AX 20	20 pz
FC1740521	NucleoBond AX 100	20 pz
FC1740525	NucleoBond AX 2000	10 pz
FC1740531	NucleoBond AX 500	10 pz
FC1740534	NucleoBond AX 10000	1 conf. da 5 pz
FC1740576	Kit NucleoBond PC2000	5 prep.
FC1740593	Kit NucleoBond PC10000	5 prep.



Accessories

Code	Description	Packaging
FC1740562	NUCLEOBOND RACK SMALL	1 pz
FC1740563	NUCLEOBOND RACK LARGE	1 pz
FC1740601	NUCLEOBOND BUFFER SET I	1 set



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

