

**NucleoBond PC EF/Prep****Macherey Nagel**

NucleoBond PC EF/Prep - per la purificazione di DNA plasmidico esente da endotossine - tecnologia: cromatografia a scambio anionico - formato: colonna - step esclusivo per la rimozione delle endotossine: non è necessaria l'incubazione refrigerata - chiarificazione del lisato mediante filtri inclusi - dimensioni del vettore 300 kbp - A260/A280 1,8/1,95 - livello endotossinico 0,1 EU/μg DNA - comprende: colonne, tamponi, RNasi A, filtri.



Code	Description	Packaging
FC1740548	Kit NucleoBond PC10000 EF	5 prep.
FC1740549	Kit NucleoBond PC2000 EF	5 prep.
FC1740550	Kit NucleoBond PC500 EF	10 prep.
FC1740594	Kit NucleoBond PC Prep 100	100 prep.

Accessories

Code	Description	Packaging
FC1740516	NUCLEOBOND BUFFER S1	500 ml
FC1740517	NUCLEOBOND BUFFER S2	500 ml
FC1740518	NUCLEOBOND BUFFER S3	500 ml
FC1740527	NUCLEOBOND BUFFER N2	500 ml
FC1740528	NUCLEOBOND BUFFER N3	1000 ml
FC1740529	NUCLEOBOND BUFFER N5	500 ml



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

