

**Protino Ni-IDA Packed Columns****Macherey Nagel**

Protino® Ni-IDA è la scelta migliore per ottenere un'alta resa in proteine His-tag.

- Gruppo chelante IDA (acido iminodiacetico), 3 siti leganti
- Alta resa/recupero anche con campioni diluiti
- Alta concentrazione proteica
- Alta specificità di legame, minore legame aspecifico con proteine contaminanti rispetto alle comuni matrici IMAC
- Resina secca: stoccaggio a temperatura ambiente
- Pronto all'uso
- Tamponi inclusi



Code	Description	Packaging
FC145150E	745150.10 PROTINO NI-IDA PACKED COL	10 prep
FC145150N	745150.50 PROTINO NI-IDA 150PACK 50	250 prep
FC1451600	PROTINO NI-IDA 100 PACK COL 5 PREP	5 prep
FC145160N	745160.50 PROT NI-IDA PACK COL 50PP	50 prep
FC1451700	745170.5 PROT NI-IDA 2000 PACK CO 5	5 prep
FC145170H	745160.25 PROT NI-IDA 2000 25 PREP	25 prep



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

