

**Protino Ni-NTA Columns 1 mL/5 mL****Macherey Nagel**

Protino® Ni-NTA, purificazione di proteine His-tag.

- Gruppo chelante NTA (acido nitrilotriacetico), 2 siti leganti
- Alta capacità ed alta affinità
- non sono necessari un'ottimizzazione e cambio di protocollo per sostituire i prodotti Ni-NTA attualmente utilizzati
- Purificazione in condizioni native e denaturanti
- Disponibile per piccole proteine, grandi complessi proteici, proteine con basso grado di espressione: uso universale
- Alto rapporto qualità/prezzo
- Pronto all'uso



Code	Description	Packaging
FC145410A	PROTINONI NTA COLUMN (5X1ML)	5X1 ml
FC145410B	PROTINONI NTA COLUMN (1X5ML)	1X5 ml
FC145410C	PROTINONI NTA COLUMN (5X5ML)	5X5 ml

Accessories

Code	Description	Packaging
FC1740681	NUCLEOVAC 96 VACUUM MANIFOLD	1pz
FC1745250	PROTINO COLUMNS 14 ml	10 colonne
FC1745255	PROTINO COLUMNS 35 ml	10 colonne



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

