



Morsetti



Code	Description	Packaging
LLG06073330	Giunto sferico cup polished KS19/9	1 pz.
LLG06073331	Giunto sferico cup polished KS29/15	1 pz.
LLG06202224	Giunto sferico cup inche KS51/30	1 pz.
LLG06202487	Giunto sferico cup polished KS64/40	1 pz.
LLG06243156	Joint clip, POM, per S 19	1 pz.
LLG06262147	Ruota dentata	1 pz.
LLG07082091	Giunto sferico cup inche KS35/25	1 pz.
LLG07082092	Giunto sferico ball inche KS35/25	1 pz.
LLG07611432	Giunto sferico cup polished KS41/25	1 pz.
LLG07611433	Giunto sferico ball polished KS19/9	1 pz.
LLG07611434	Giunto sferico ball polished KS35/20	1 pz.
LLG07611435	Giunto sferico ball polished KS41/25	1 pz.