



Pallone fondo tondo



Code	Description	Packaging
LLG07021822	Pallone fondo piatto, 250 mL, NS 24/29, DURAN	1 pz.
LLG07021824	Pallone fondo piatto, 500 mL, NS 24/29, DURAN	1 pz.
LLG07021829	Pallone fondo tondo, 100 mL, 2 colli, NS24/28	1 pz.
LLG07021831	Pallone fondo tondo, 250 mL, 2 colli, NS24/29	1 pz.
LLG07021833	Pallone fondo tondo, 500 mL, 2 colli, NS24/30	1 pz.
LLG07021835	Pallone fondo tondo, 1000 mL, 2 colli, NS24/31	1 pz.
LLG07021837	Pallone fondo tondo, 2000 mL, 2 colli, NS24/32	1 pz.
LLG07021849	Pallone fondo tondo, 250 mL, 2 colli, NS24/33	1 pz.
LLG07021851	Pallone fondo tondo, 500 mL, 3 colli, NS24/29	1 pz.
LLG07021853	Pallone fondo tondo, 1000 mL, 3 colli, NS24/30	1 pz.

