



Accessori per Titronic® 300 e TitroLine® 5000

SI Analytics

Accessori per Titronic® 300 e TitroLine® 5000.



Code	Description	Packaging
LLG09863188	Agitatore magnetico TM 50;	1 pz.
LLG09863190	Bottiglia ambrata, 1000ml, tipo TZ 3803;	1 pz.
LLG09863189	Connessione USB-HUB tipo TZ 3830;	1 pz.



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

