



**Accessori per Titronic® 300 e TitroLine® 5000**

**SI Analytics**

Accessori per Titronic® 300 e TitroLine® 5000.



| Code        | Description                              | Packaging |
|-------------|------------------------------------------|-----------|
| LLG09863188 | Agitatore magnetico TM 50;               | 1 pz.     |
| LLG09863190 | Bottiglia ambrata, 1000ml, tipo TZ 3803; | 1 pz.     |
| LLG09863189 | Connessione USB-HUB tipo TZ 3830;        | 1 pz.     |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

