



Setti per tappi a vite ND9

LLG Labware

Setti per tappi a vite ND9.



Code	Description	Packaging
LLG04008213	Setto N9 gomma rossa/FEP incolore durezza 40°/A spessore 1mm	100 pz
LLG04008211	Setto N9 PTFE vergine, bianco durezza 53°/D spessore 0,25 mm	100 pz
LLG04008212	Setto N9 silicone bianco/PTFE rosso durezza 40°/A spessore 1 mm	100 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

