

**Provette per centrifuga con fondo tondo DURAN®****DURAN Group**

Provette per centrifuga con fondo tondo DURAN®. DIN 58970 Parte 2. Vetro tipo 1 / vetro neutro, in conformità con USP, EP e JP. Autoclavabili. RCF max 4000 x g.



Code	Description	Packaging
LLG09315011	Provette per centrifuga, fondo tondo, DURAN, 12 ml, Ø 16 mm, h 100 mm	1 pz
LLG09315036	Provette per centrifuga, fondo tondo, DURAN, 250* ml, Ø 56 mm, h 147 mm	1 pz
LLG09315014	Provette per centrifuga, fondo tondo, DURAN, 25 ml, Ø 24 mm, h 100 mm	1 pz
LLG09315017	Provette per centrifuga, fondo tondo, DURAN, 50 ml, Ø 34 mm, h 100 mm	1 pz
LLG09315010	Provette per centrifuga, fondo tondo, DURAN, 6 ml, Ø 12 mm, h 100 mm	1 pz
LLG09315024	Provette per centrifuga, fondo tondo, DURAN, 80* ml, Ø 40 mm, h 115 mm	1 pz
LLG09315026	Provette per centrifuga, fondo tondo, DURAN, 80 ml, Ø 44 mm, h 100 mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

