



Mestoli per latte

Gerber

Mestoli per latte. In ottone nichelato, con valvola, sterilizzabili.



Code	Description	Packaging
LLG06266491	Mestolo per latte, cap. 1 ml	1 pz.
LLG06266495	Mestolo per latte, cap. 10 ml	1 pz.
LLG06266499	Mestolo per latte, cap. 100 ml	1 pz.
LLG06266492	Mestolo per latte, cap. 2 ml	1 pz.
LLG06266496	Mestolo per latte, cap. 20 ml	1 pz.
LLG06266493	Mestolo per latte, cap. 4 ml	1 pz.
LLG06266497	Mestolo per latte, cap. 40 ml	1 pz.
LLG06266494	Mestolo per latte, cap. 5 ml	1 pz.
LLG06266498	Mestolo per latte, cap. 50 ml	1 pz.

