

**Sonde per omogenizzatori ad Ultrasuoni Sonopuls, titanio****Bandelin**

Sonde per omogenizzatori ad Ultrasuoni Sonopuls, titanio. Sonde in titanio. Per collegamento a sonde pluristadio / sonde ausiliarie.



| Code        | Description                                                                 | Packaging |
|-------------|-----------------------------------------------------------------------------|-----------|
| LLG06274618 | Sonde per Sonopuls HD 4400, Tipo TS 413, Per SH 400 G, 100-750 ml, Ø 13mm   | 1 pz      |
| LLG09650169 | Sonde per Sonopuls HD 4400, Tipo TS 416, Per SH 400 G, 250-1000 ml, Ø 16mm  | 1 pz      |
| LLG07647594 | Sonde per SH 100 G SH 200 G, Per 20-200 20-900 ml, Per HD 4100/4200, Ø13mm  | 1 pz      |
| LLG09650239 | Sonda tipo VS 70 T, Ø 13 mm, puntale lungo, per HD 2070/2200 HD 3100/3200   | 1 pz.     |
| LLG09650247 | Sonda tipo TT 13 FZ, Ø 13mm, piastra titanio c/perforaz., per HD 2200/3200  | 1 pz.     |
| LLG09650250 | Sonda tipo TT 13, Ø13mm, piastra in titanio, per HD 2070/2200 HD 3100/3200  | 1 pz.     |
| LLG09650243 | Sonda tipo VS 190 T, Ø 19,0 mm, puntale lungo, per HD 2200/3200HD 3400      | 1 pz.     |
| LLG06251521 | Sonda tipo MS 2.0, Ø 2,0mm, punta Micro, per 0,25-20ml, per Sonpuls mini20  | 1 pz.     |
| LLG09650157 | Sonda in titanio TS 102, Ø 2,5mm, per SH 100 G                              | 1 pz.     |
| LLG09650244 | Sonda tipo MS 72, Ø 2mm, Micropunta, per HD 2070/2200 HD 3100/3200          | 1 pz.     |
| LLG09650245 | Sonda tipo MS 73, Ø 3mm, Micropunta, per HD 2070/2200 HD 3100/3200          | 1 pz.     |
| LLG09650159 | Sonda in titanio TS 104, Ø 4mm, per SH 100 G                                | 1 pz.     |
| LLG09650161 | Sonde x SH 100 G SH 200 G, x 5-75 10-100 10-350ml, x HD 4050/4100/4200 Ø6mm | 1 pz      |
| LLG09650246 | Sonda tipo KE 76, Ø 6 mm, punta conica, per HD 2070/2200 HD 3100/3200       | 1 pz.     |
| LLG09650240 | Sonda tipo VS 200 T, Ø 25,0 mm, puntale lungo, per HD 2200/3200 HD 3400     | 1 pz.     |
| LLG09650249 | Sonda tipo TT 25, Ø 25mm, piastra in titanio, per HD 2200/3200              | 1 pz.     |
| LLG09650165 | Sonda tipo MS 1.5, Ø 1,5mm, punta Micro, per 0,1-10 ml, per Sonpuls mini20  | 1 pz.     |
| LLG09650166 | Sonda tipo MS 2.5, Ø 2,5 mm, punta Micro, per 0,5-25 ml, per Sonpuls mini20 | 1 pz.     |
| LLG09650158 | Sonda in titanio TS 103, Ø 3mm, per SH 100 G                                | 1 pz.     |





**CARLO ERBA Reagents S.r.l.**  
VIA R. MERENDI, 22  
20007 CORNAREDO (MI)  
TEL. +39 02 93 99 190  
FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
a Certified Quality  
Management System*

