

**Sonda esterna a 2 canali per USB data Logger EBI 310****ebro**

Sonda esterna a 2 canali per USB data Logger EBI 310. Adattatore per EBI 310 per misurare due temperature parallele. Si prega di ordinare 2 sensori separatamente.



Code	Description	Packaging
LLG06264939	Sonda per temperatura TPX 310-P1	1 pz.
LLG06264938	Sonda esterna per temperatura TPX 310	1 pz.



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

