



Flask on slide Lab-Tek

Thermo Scientific

Flask on slide Nunc# Lab-Tek# - sterile - le cellule crescono direttamente sul vetrino da microscopio - ideale per l'autoradiografia o l'immunofluorescenza di singole cellule - VERSIONE Slide Flask: vetrino in polistirolo trattato per colture cellulari saldato mediante ultrasuoni a una flask dim. 18x50 mm, tappo nero, la struttura superiore si separa facilmente dal vetrino dopo la coltura liberandolo per l'analisi al microscopio, certificata Nunclon# - VERSIONE Flaskerette chamber: vetrino in vetro con guarnizione in silicone e flask dim. 20x50 mm, può essere sigillata per impedire il passaggio di aria, ideale per la preparazione di cellule da analizzare al microscopio elettronico - marcato CE.



Code	Description	Packaging
FA9170920	flask on slide SlideFlask - in PS trattato - 9,0 cm2	50 pezzi
FA9177453	flask on slide Flaskette - in vetro - 10 cm2	96 pezzi

Accessories

Code	Description	Packaging
FA9171862	coprivetrino per FA9170920 - 18x50 mm	1000 pezzi



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

