



Lastra per TLC, strato gel di silice SIL HD

Macherey Nagel

Lastra per TLC, strato gel di silice SIL HD.



Code	Description	Packaging
MNL809210	TLC plates SIL HD, 10x10 cm	25 pz
MNL809212	TLC plates SIL HD, 10x20 cm	50 pz
MNL809213	TLC plates SIL HD, 20x20 cm	25 pz
MNL809217	TLC plates SIL HD, 5x10 cm	50 pz
MNL809220	TLC plates SIL HD UV254, 10x10 cm	25 pz
MNL809222	TLC plates SIL HD UV254, 10x20 cm	50 pz
MNL809223	Lastra TLC, gel di silice SIL HD, UV254, 20x20 cm	25 lastre
MNL809227	TLC plates SIL HD UV254, 5x10 cm	50 pz

