


**Burette di ricambio per burette automatiche di Schilling SILBERBRAND****BRAND**

Burette di ricambio per burette automatiche di Schilling, SILBERBRAND. Disponibili in boro 3.3 o AR-GLAS®.



Code	Description	Packaging
FB423723	Buretta ricambio burette Schilling, 25 ml, vetro ambrato, 390 mm, AR-GLAS	1 pz.



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

