

**Starter-Kit F1-ClipTip****Thermo Scientific**

Starter-Kit F1-ClipTip - da utilizzare con puntali Clip-Tip - il KIT 1 comprende: micropipetta monocanale F1-ClipTip 1-10 μ L, micropipetta monocanale F1-ClipTip 10-100 μ L, micropipetta monocanale F1-ClipTip 100-1000 μ L, stand, puntali ClipTip 20 μ L in rack, puntali ClipTip 200 μ L in rack, puntali ClipTip 1000 μ L in rack, vaschetta da 25 mL, vaschetta a 2 scomparti da 25 mL, vaschetta a 2 scomparti da 100 mL - il KIT 2 comprende: micropipetta monocanale F1-ClipTip 2-20 μ L, micropipetta monocanale F1-ClipTip 20-200 μ L, micropipetta monocanale F1-ClipTip 100-1000 μ L, micropipetta a 8 canali F1-ClipTip 30-300 μ L, stand, puntali ClipTip 20 μ L in rack, puntali ClipTip 200 μ L in rack, puntali ClipTip 300 μ L in rack, puntali ClipTip 1000 μ L in rack, vaschetta da 25 mL, vaschetta a 2 scomparti da 25 mL, vaschetta a 2 scomparti da 100 mL.



Code	Description	Packaging
FL7701140	starter-kit GLP kit 1	1 pezzo
FL7701150	starter-kit GLP kit 2	1 pezzo



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

