



Testate per pompe a pistone in ceramica MCP-CPF, senza valvole

Ismatec

Testate per pompe a pistone in ceramica MCP-CPF, senza valvole.



Code	Description	Packaging
LLG09828430	Testate QP.QO.SSY pist. acciaio SS316 3,2-80µl 0,04-144ml/min Øtubo 6,4mm	1 pz.
LLG09828424	Testate QP.Q1.CSC-WT** pist. ceramica 12,8-320µl 0,13-576ml/min Øtubo 9,5mm	1 pz.
LLG09828425	Testate QP.Q2.CSC-WT**pist. ceramica 28,8-720µl 0,29-1300ml/min Øtubo 9,5mm	1 pz.
LLG09828432	Testate QP.Q2.CSC pist. ceramica 28,8-720µl 0,29-1300ml/min Øtubo 9,5mm	1 pz.
LLG09828433	Testate QP.Q3.CKC pist.ceramica 51,2-1280 µl 0,51-2300 ml/min Øtubo 12,7 mm	1 pz.



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

