

**Imbuti filtranti, Buchner, DURAN®****DURAN Group**

Imbuti filtranti, Buchner, DURAN®. Vetro DURAN®. Vetro tipo 1 / vetro neutro, in conformità con USP, EP e JP. Autoclavabile.



Code	Description	Packaging
LLG09052341	Imbuti Buchner, 1000 ml, porosità 1, Ø 120 mm, Ø stelo 22 mm, h 270 mm	1 pz
LLG09052342	Imbuti Buchner, 1000 ml, porosità 2, Ø 120 mm, Ø stelo 22 mm, h 270 mm	1 pz
LLG09052343	Imbuti Buchner, 1000 ml, porosità 3, Ø 120 mm, Ø stelo 22 mm, h 270 mm	1 pz
LLG09052344	Imbuti Buchner, 1000 ml, porosità 4, Ø 120 mm, Ø stelo 22 mm, h 270 mm	1 pz
LLG09052345	Imbuti Buchner, 1000 ml porosità 5 Ø 120 mm Ø stelo 22 mm h 320 mm	1 pz
LLG09052321	Imbuti Buchner, 125 ml, porosità 1, Ø 60 mm, Ø stelo 10 mm, h 140 mm	1 pz
LLG09052322	Imbuti Buchner, 125 ml, porosità 2, Ø 60 mm, Ø stelo 10 mm, h 140 mm	1 pz
LLG09052323	Imbuti Buchner, 125 ml, porosità 3, Ø 60 mm, Ø stelo 10 mm, h 140 mm	1 pz
LLG09052324	Imbuti Buchner, 125 ml, porosità 4, Ø 60 mm, Ø stelo 10 mm, h 140 mm	1 pz
LLG09052325	Imbuti Buchner, 125 ml porosità 5 Ø 60 mm Ø stelo 10 mm h 140 mm	1 pz
LLG09052361	Imbuti Buchner, 4000 ml, porosità 1, Ø 175 mm, Ø stelo 30 mm, h 425 mm	1 pz
LLG09052362	Imbuti Buchner, 4000 ml, porosità 2, Ø 175 mm, Ø stelo 30 mm, h 425 mm	1 pz
LLG09052363	Imbuti Buchner, 4000 ml, porosità 3, Ø 175 mm, Ø stelo 30 mm, h 425 mm	1 pz
LLG09052364	Imbuti Buchner, 4000 ml, porosità 4, Ø 175 mm, Ø stelo 30 mm, h 425 mm	1 pz
LLG09052301	Imbuti Buchner, 50 ml, porosità 1, Ø 35 mm, Ø stelo 10 mm, h 130 mm	1 pz
LLG09052302	Imbuti Buchner, 50 ml, porosità 2, Ø 35 mm, Ø stelo 10 mm, h 130 mm	1 pz
LLG09052303	Imbuti Buchner, 50 ml, porosità 3, Ø 35 mm, Ø stelo 10 mm, h 130 mm	1 pz
LLG09052304	Imbuti Buchner, 50 ml, porosità 4, Ø 35 mm, Ø stelo 10 mm, h 130 mm	1 pz
LLG09052305	Imbuti Buchner, 50 ml porosità 5 Ø 35 mm Ø stelo 10 mm h 130 mm	1 pz
LLG09052331	Imbuti Buchner, 500 ml, porosità 1, Ø 95 mm, Ø stelo 22 mm, h 240 mm	1 pz



Code	Description	Packaging
LLG09052332	Imbuti Buchner, 500 ml, porosità 2, Ø 95 mm, Ø stelo 22 mm, h 240 mm	1 pz
LLG09052333	Imbuti Buchner, 500 ml, porosità 3, Ø 95 mm, Ø stelo 22 mm, h 240 mm	1 pz
LLG09052334	Imbuti Buchner, 500 ml, porosità 4, Ø 95 mm, Ø stelo 22 mm, h 240 mm	1 pz
LLG09052335	Imbuti Buchner, 500 ml porosità 5 Ø 95 mm Ø stelo 22 mm h 200 mm	1 pz
LLG09052311	Imbuti Buchner, 75 ml, porosità 1, Ø 45 mm, Ø stelo 10 mm, h 132 mm	1 pz
LLG09052312	Imbuti Buchner, 75 ml, porosità 2, Ø 45 mm, Ø stelo 10 mm, h 132 mm	1 pz
LLG09052313	Imbuti Buchner, 75 ml, porosità 3, Ø 45 mm, Ø stelo 10 mm, h 132 mm	1 pz
LLG09052314	Imbuti Buchner, 75 ml, porosità 4, Ø 45 mm, Ø stelo 10 mm, h 132 mm	1 pz
LLG09052315	Imbuti Buchner, 75 ml porosità 5 Ø 45 mm Ø stelo 10 mm h 130 mm	1 pz

