

**Filtri/setacci cellulari, cell strainer, in nylon, sterili****JetBiofil**

Filtri cellulari, Nylon, sterili. I filtri cellulari sono particolarmente adatti per la generazione di sospensioni omogenee unicellulari da tessuti o cellule primarie di lisati.

- Disponibili in 3 misure mesh 40 μm , 70 μm e 100 μm
- Codice colore per la facile determinazione della misura delle mesh.
- Studiati per adattarsi alle provette Jet Biofil coniche da 50 ml.
- Confezionati individualmente.
- La linguetta superiore permette la manipolazione asettica mediante pinzette.
- Mesh in nylon altamente resistente con pori uniformemente distanziati.
- Apirogeni.
- Esenti da DNasi- e RNasi.
- Sterilizzati con raggi gamma.
- Proteggono la strumentazione per citometria di flusso e cell sorting mediante l'affidabile rimozione di grumi (clump) e pezzi di cellule.



Code	Description	Packaging
CSS013040	Cell Strainer, in nylon, blu, 40 µm, sterile	50 pz
CSS013070	Cell Strainer, in nylon, bianco, 70 µm, sterile	50 pz
CSS013100	Cell Strainer, in nylon, giallo, 100 µm, sterile	50 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

