



Flaconi da centrifuga bocca larga tipo 3141, in PP-copolimero

Flaconi da centrifuga bocca larga tipo 3141, in PP-copolimero. In PP copolimero con tappo a vite di tenuta. Autoclavabili.



Code	Description	Packaging
LLG09315763	Flaconi da centrifuga 3141, Ø 60,7 mm, h 133,4 mm, max. RCF 27500 x g	1 pz.
LLG09315764	Flaconi da centrifuga 3141, Ø 69,5 mm, h 160,0 mm, max. RCF 13700 x g	1 pz.



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

