

**Accessori per sonopuls HD 2200, in vetro****Bandelin**

Accessori per sonopuls HD 2200, in vetro. Celle per trattamento a ultrasuoni diretto. In vetro.



Code	Description	Packaging
LLG09650261	Cella Rosette per sonotrodi Ø max.13mm, 110 ml, per HD 2070/2200/3100/3200	1 pz.
LLG09650235	Cella con camicia raffreddamento, Ø13mm, 15ml, per HD 2070/2200/3100/3200	1 pz.
LLG09650236	Cella flusso continuo c/camicia raffr. Ø13mm 15ml, x HD2070/2200/3100/3200	1 pz.
LLG09650223	Cella con camicia raffreddamento, Ø25mm, 70ml, per HD 2200/3200	1 pz.
LLG09650267	Cella Rosette per sonotrodi Ø max.13 mm, 25 ml, per HD 2070/2200/3100/3200	1 pz.
LLG09650260	Cella Rosette per sonotrodi Ø max. 25 mm, 390 ml, per HD 2200/3200/3400	1 pz.
LLG09650268	Cella Rosette per sonotrodi Ø max.13 mm, 40 ml, per HD 2070/2200/3100/3200	1 pz.
LLG09650264	Cella Rosette per sonotrodi Ø max. 25 mm, 660 ml, per HD 2200/3200/3400	1 pz.
LLG09650224	Cella a flusso continuo c/camicia raffredd., Ø25mm, 70ml, per HD 2200/3200	1 pz.