



Adattatori per centrifughe Hermle

HERMLE Labortechnik GmbH

Adattatori per centrifughe Hermle



Code	Description	Packaging
LLG04682541	Adattatore per 1 x 6 ... 10 ml Vacutainer® e provette per prelievo ematico da 9 ml, Ø 17mm, max. L 107mm per 221.54 V01 / 221.54 V04 / 221.54 V20 / 221.57 V02 / 221.19 V20	1 pz
LLG06225443	Adattatore per 1 provetta da 50 ml tubes (conica), per 220.78 V21 / 221.18 V20 / 221.18 V21 / 221.20 V20, conf. 2 pz	2 pz
LLG06225444	Adattatore per 1 provetta da 15 ml (conica), per 220.78 V21 / 221.18 V20 / 221.18 V21 / 221.20 V20, conf. 2 pz	2 pz
LLG06237852	Adattatore per 1 provetta x 15 ml, per 221.22 V20 / 221.22 V21 / 221.55 V01 / 221.55 V04 / 221.55 V20 , conf. 2 pz	1 pz
LLG06238600	Adattatore per 1 provetta x 50 ml, per 220.78 V21 / 221.18 V20 / 221.18 V21 / 221.20 V20, conf. 2 pz	1 pz
LLG06239997	Adattatore per 1 provetta x 15 ml, per 220.78 V21 / 221.18 V20 / 221.18 V21 / 221.20 V20, conf. 2 pz	1 pz
LLG06239998	Adattatore per 1 provetta da 30 ml, Ø 25 mm, lungh. 90/103 mm, per 220.78 V21 / 221.18 V20 / 221.18 V21 / 221.20 V20, conf. 2 pz	1 pz
LLG06239999	Adattatore per 2 provette da 15 ml, Ø 17 mm, lungh. 100 mm, per 220.78 V21 / 221.18 V20 / 221.18 V21 / 221.20 V20, conf. 2 pz	1 pz
LLG06240000	Adattatore per 1 provetta da 16 ml, Ø 18mm, lungh. 103/110 mm, per 220.78 V21 / 221.18 V20 / 221.18 V21 / 221.20 V20, conf. 2 pz	1 pz
LLG06258502	Set adattatore per 1.5 / 2.0 ml contenitori di reazione (per RCF > 16000 x g), per 221.23 V20, conf. 6 pz	1 pz
LLG06261823	Set adattatori per contenitori di reazione 1.5/2.0 ml	1 pz
LLG06290008	Adattatore per 1 x 16 ml dia. 18 mm, lungh. min./max. 100/105 mm, conf. 2	1 pz
LLG07980344	Adattatore per 1 provetta x 30 ml, per 221.55 V01 / 221.55 V04 / 221.55 V20 / 221.22 V20 / 221.22 V21, conf. 2 pz	1 pz
LLG07980345	Adattatore per 1 provetta x 15 ml (conica), per 221.22 V20 / 221.22 V21, conf. 2 pz	2 pz





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

