



Provetta per centrifuga in vetro DURAN®

Glaswarenfabrik Karl Hecht

Centrifuge tubes, conical, with rim, borosilicate glass 3.3

- Fused white graduation
- Mark at 15 ml.



Code	Description	Packaging
LLG06244140	Provetta da centrifuga 112x17mm, in vetro DURAN®, conica, con bordo, graduata, cap. 10 ml, grad. 0.1 ml	1 pz
LLG06281060	Provetta per centrifuga in vetro DURAN®, 15 ml, Ø 17 mm, H 112 mm, fondo conico, con bordo, graduata fino a 15 ml, grad. 0,1 ml	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

