

**Piastra PCR FrameStar Break-A-Way****AZENTA**

Piastra PCR FrameStar Break-A-Way.

- 96 pozzetti. Semi-skirted;
- Pre-intaglio verticale per permettere di dividerla in strip da 8 o parti di piastra;
- Angolo tagliato $\geq 12^\circ$;
- Profilo standard;
- Pozzetti in polipropilene ≥ 0.2 ml;
- Cornice in policarbonato;
- Volume di lavoro: $\geq 200 \mu\text{l}$;
- Capacit  totale del pozzetto: $300 \mu\text{l}$;
- Esente da DNasi, RNasi, DNA genomico umano, pirogeni e endotossine;
- Adatta alla maggior parte dei termociclatori;
- Compatibile con le pipette multicanali standard.

Pozzetti

Pareti ultra-sottili e uniformi per risultati ottimali in PCR e real-time (RT-qPCR).

Cornice

Linguette terminali per facilitare la manipolazione e l'etichettatura. Griglia alfanumerica per la facile identificazione del campione.

Caratteristiche

- Lungh. piastra: 125.11 ± 0.25 mm;
- Largh. piastra: 83.22 ± 0.25 mm;
- Altezza piastra: 20.80 ± 0.25 mm;
- Profondit  pozzetto: 20.30 ± 0.15 mm;
- Diam. pozzetto: 5.46 ± 0.10 mm;
- Capacit  totale del pozzetto: 0.3 ml ($300 \mu\text{l}$);
- Volume di lavoro ≥ 0.2 ml ($200 \mu\text{l}$).



Code	Description	Packaging
FE94TI1000B	Break-A-Way 96 pozz. semi-skirt, scanalature verticali, frame PC blu, A12	50 plates
FE94TI1000C	Break-A-Way 96 pozz. semi-skirt, scanalature verticali, frame PC trasp., A1	50 plates
FE94TI1000G	Break-A-Way 96 pozz. semi-skirt, scanalature verticali, frame PC verde, A12	50 plates
FE94TI1000P	Break-A-Way 96 pozz. semi-skirt, scanalature verticali, frame PC viola, A12	50 plates
FE94TI1000R	Break-A-Way 96 pozz. semi-skirt, scanalature verticali, frame PC rosso, A12	50 plates
FE94TI1000W	Break-A-Way 96 pozz. semi-skirt, scanalature verticali, frame PC bianco, A1	50 plates
FE94TI1000X	Break-A-Way 96 pozz. semi-skirt, scanalature verticali, frame PC nero, A12	50 plates
FE94TI1001	Break-A-Way 96 pozz. semi-skirt, scanalature verticali, frame PC nero, A12	50 plates
FE94TI1200B	Break-A-Way 96 semi-skirt LP, scanalature verticali, frame PC blu, A12	50 plates
FE94TI1200C	Break-A-Way 96 semi-skirt LP, scanalature verticali, frame PC trasp., A12	50 plates
FE94TI1200G	Break-A-Way 96 semi-skirt LP, scanalature verticali, frame PC verde, A12	50 plates
FE94TI1200P	Break-A-Way 96 semi-skirt LP, scanalature verticali, frame PC viola, A12	50 plates
FE94TI1200R	Break-A-Way 96 semi-skirt LP, scanalature verticali, frame PC rosso, A12	50 plates
FE94TI1200W	Break-A-Way 96 semi-skirt LP, scanalature verticali, frame PC bianco, A12	50 plates
FE94TI1200X	Break-A-Way 96 semi-skirt LP, scanalature verticali, frame PC nero, A12	50 plates
FE94TI1201	Break-A-Way 96 bianchi semi-skirt LP scanalature verticali, frame nero, A12	50 plates



CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

