

**Strip per PCR Vari-Strip, 8 pozzetti****AZENTA**

Strip per PCR Vari-Strip[®], 8 pozzetto, Low Profile.

- Pozzetti in polipropilene da 0,1 mL.
- Volume di lavoro $< 100 \mu\text{l}$.
- Capacit  totale del pozzetto: $200 \mu\text{l}$.
- Strip utilizzabili da sole o in combinazione con cornici Vari-Plate[®] per una piastre; facile manipolazione.
- Esenti da DNasi, RNasi, DNA genomico umano, endotossine, pirogeni.



Code	Description	Packaging
FE94TI0753	Vari-Strip 8 pozz. Low Profile, PP bianco	120 tube strips
FE94TI0753C	Vari-Strip 8 pzz. Low Profile PP trasparente	120 tube strips
FE94TI0754	Vari-Strip 8 pzz. Low Profile PP bianco, con strip tappi piatti	120 tube strips + cap strips
FE94TI0754C	Vari-Strip 8 pzz. Low Profile PP trasparente, con strip tappi piatti	120 tube strips + cap strips



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

