



Bacinelle per laboratorio in PS

Vassoi e vaschette, PS. Adatti per usi alimentari ed usi generali di laboratorio.
Polistirene antiurto.



Code	Description	Packaging
LLG06233344	Vassoio, dim. 201x151x21mm, dim. interne 170x110x10 mm	1 pz.
LLG06228312	Vassoi e vaschette, PS, Largh. 202 mm, Prof. 151 mm, Altezza 81 mm, Dim. interne (L x P x H) 140 x 90 x 70 mm	1 pz
LLG06224986	*Tray for suitable foodstuffs	1 pz
LLG06225055	Vassoi e vaschette, PS, Largh. 353 mm, Prof. 253 mm, Altezza 41 mm, Dim. interne (L x P x H) 310 x 210 x 30 mm	1 pz.
LLG06206152	Vassoio, Largh. 254 mm, Profondità 254 mm, Altezza 21 mm	1 pz.
LLG06225054	Vassoi e vaschette, PS, Largh. 303 mm, Prof. 151 mm, Altezza 42 mm, Dim. interne (L x P x H) 260 x 110 x 30 mm	1 pz.
LLG06224951	Vassoio in PS, Largh. 409 mm, Prof. 300 mm, Altezza 42 mm, Dim. interne (L x P x H) 360 x 250 x 30 mm	1 pz.
LLG06227810	Vassoi e vaschette, PS, Largh. 352 mm, Prof. 252 mm, Altezza 81 mm, Dim. interne (L x P x H) 290 x 190 x 70 mm	1 pz.
LLG06205063	Vassoi e vaschette, PS, Largh. 408 mm, Prof. 300 mm, Altezza 21 mm, Dim. interne (L x P x H) 380 x 270 x 10 mm	1 pz
LLG06227811	Vassoi e vaschette, PS, Largh. 408 mm, Prof. 299 mm, Altezza 81 mm, Dim. interne (L x P x H) 340 x 240 x 70 mm	1 pz





CARLO ERBA Reagents S.r.l.
VIA R. MERENDI, 22
20007 CORNAREDO (MI)
TEL. +39 02 93 99 190
FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
a Certified Quality
Management System*

