

**Termometri a resistenza ScienceLine**

Xylem Analytics Germany (SI)

Termometri a resistenza ScienceLine

Range temperatura da -30 a 135 °C, con sensore Pt 1000, cavo di collegamento Ø 2 x 4 mm (W 5790 PP con Ø 2 x 2 mm), cavo da 1 mt (W 2180-KOAX con KOAX testa).



Code	Description	Packaging
LLG09041223	Termometri a resistenza ScienceLine, Tipo W 5790 NN, Lungh. 120 mm, Ø 4 mm, Materiale stelo Acciaio inox	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

