



Tubo per laboratorio Isoflex

Tubo per laboratorio Isoflex. Tubo in PVC atossico, autoclavabile, trasparente, anti-invecchiamento e anti-ossidazione. Imballato in sacchetti di PE. Il tubo è marcato ad ogni metro. Buona resistenza ad acidi e basi forti, scarsa resistenza ai solventi organici. Durometro 77 Shore A. Confezionato in rotolo da 20 m.



Code	Description	Packaging
LLG09205147	Tubo in PVC per laboratorio, Øint. 7 mm, Ø est. 11,0 mm, spess. 2,0 mm	20 pz
LLG09205148	Tubo in PVC per laboratorio, Øint. 8 mm, Ø est. 12,0 mm, spess. 2,0 mm	20 pz
LLG09205149	Tubo in PVC per laboratorio, Øint. 10 mm, Ø est. 14,0 mm, spess. 2,0 mm	20 pz
LLG09205150	Tubo in PVC per laboratorio, Øint. 12 mm, Ø est. 17,0 mm, spess. 2,5 mm	20 pz
LLG09205142	Tubo in PVC per laboratorio, Øint. 1 mm, Ø est. 2,4 mm, spess. 0,8 mm	20 pz
LLG09205143	Tubo in PVC per laboratorio, Øint. 2 mm, Ø est. 4,0 mm, spess. 1,0 mm	20 pz
LLG09205144	Tubo in PVC per laboratorio, Øint. 4 mm, Ø est. 6,0 mm, spess. 1,0 mm	20 pz
LLG09205145	Tubo in PVC per laboratorio, Øint. 5 mm, Ø est. 8,0 mm, spess. 1,5 mm	20 pz
LLG09205146	Tubo in PVC per laboratorio, Øint. 6 mm, Ø est. 9,0 mm, spess. 1,5 mm	20 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

