



Accessori per omogenizzatore BagSystem®

interscience

Accessori per omogenizzatore BagSystem®.



Code	Description	Packaging
LLG09570088	Clip per chiudere i sacchetti BagClip® 80/100	50 pz.
LLG09570089	Clip per chiudere i sacchetti BagClip® 400	50 pz.
LLG09570090	Clip per chiudere i sacchetti BagClip® 3500	50 pz.
LLG09570083	Apri-sacchetti BagOpen® 80/100	1 pz.
LLG09570084	Apri-sacchetti BagOpen® 400	1 pz.
LLG09570085	Apri-sacchetti BagOpen® 3500	1 pz.
LLG09570087	Sostegno per 10 sacchetti BagRack® 400	1 pz.
LLG09570086	Sostegno per 11 sacchetti BagRack® 80/100	1 pz.
LLG09570091	Pipetta con 3 graduazioni fisse: 0,1-0,9-1,0 ml BagPipet®	1 pz.
LLG09570092	Puntali BagTips® Regular, per BagPipet®, 19 cm, sterile	1000 pz.
LLG09570093	Puntali BagTips® Jumbo, per BagPipet®, 24 cm, sterile	1000 pz.
LLG06266940	BagOpen®, strisce adesive (set)	10 pz

