



**Test in cuvetta per analisi dell'acqua**

**HACH LANGE**

Test in cuvetta per analisi dell'acqua



| Code        | Description                                                                            | Packaging |
|-------------|----------------------------------------------------------------------------------------|-----------|
| LLG06803024 | Tensioattivi anionici, test in cuvetta, 0,05-2,0 mg/L                                  | 1 pz      |
| LLG06802583 | Test in cuvetta per COD, 100 - 2000 mg/l O <sub>2</sub> , 25 test                      | 1 pz      |
| LLG06231183 | Test in cuvetta ammonio 2,0-47,0 mg/l NH <sub>4</sub> -N, 25 test                      | 1 pz      |
| LLG06231248 | Test in cuvetta per cloruro, 1 - 70 mg/l / 70 - 1000 mg/l Cl, 24 test                  | 1 pz      |
| LLG06231242 | Test in cuvetta per azoto totale 20-100 mg/L TNb, 25 test                              | 1 pz      |
| LLG06231184 | Cuvette Tipo LCK349, per fosfato totale, 0.05 - 1.50mg/L PO <sub>4</sub> -P            | 1 pz      |
| LLG06242605 | Test in cuvetta per fosfato (orto/totale), 2,0 - 20,0 mg/l PO <sub>4</sub> -P, 25 test | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

