



Mortai in porcellana

LLG Labware

Mortai in porcellana. Mortaio privo di pestello.



Code	Description	Packaging
LLG06227289	Mortaio, 110 ml, Ø est. 83 mm, h 50 mm, lucidato	1 pz
LLG06232577	Mortaio, 110 ml, Ø est. 83 mm, h 50 mm, interno ruvido	1 pz
LLG06243865	Mortaio, 160 ml, Ø est. 90 mm, h 56 mm, DIN 12906, interno ruvido	1 pz
LLG06234591	Mortaio, 220 ml, Ø est. 105 mm, h 64 mm, lucidato	1 pz
LLG09164412	Mortaio, 220 ml, Ø est. 105 mm, h 64 mm, interno ruvido	1 pz
LLG06252383	Mortaio, 2250 ml, Ø est. 224 mm, h 127 mm, interno ruvido	1 pz
LLG09164410	Mortaio, 25 ml, Ø est. 51 mm, h 32 mm, interno ruvido	1 pz
LLG09164413	Mortaio, 400 ml, Ø est. 125 mm, h 71 mm, DIN 12906, interno ruvido	1 pz
LLG06227149	Mortaio, 70 ml, Ø est. 63 mm, h 41 mm, DIN 12906, interno ruvido	1 pz
LLG09164411	Mortaio, 75 ml, Ø est. 76 mm, h 45 mm, interno ruvido	1 pz

