

**Provette per centrifuga, coniche, PP/HDPE****Thermo Scientific**

Provette per centrifuga, coniche, PP/HDPE. Nuove e migliorate provette coniche con la più larga area di scrittura presente sul mercato, per incrementare la tracciabilità grazie allo spazio aggiunto, per inserimento di informazioni sull'etichetta. Il supporto plastico riciclabile riduce lo scarto in laboratorio. Le nuove provette coniche hanno maggior rapporto RCF, che consente un numero maggiore di applicazioni per centrifugazioni da velocità bassa a molto elevata. Sono conformi ai requisiti USP classe VI e sono non-citotossiche, apirogene e prive di RNase/DNase. Graduate, a tenuta e sterili (SAL 10-6) Materiale provette: PPMateriale tappi: HDPE



Code	Description	Packaging
FA9337846	Provette coniche, PP/HDPE, 11 ml, H 110 mm, Max. 3000xg	480 pezzi
FA9339650	Provette coniche, PP/HDPE, Sfuse, 15 ml, Ø 17 mm, H 120 mm, Max. 10500xg	500 pezzi
FA9339651	Provette coniche, PP/HDPE, In Rack, 15 ml, Ø 17 mm, H 120 mm, Max. 10500xg	500 pezzi
FA9339652	Provette coniche, PP/HDPE, Sfuse, 50 ml, Ø 30 mm, H 115 mm, Max. 17000xg	500 pezzi
FA9339653	Provette coniche, PP/HDPE, In Rack, 50 ml, Ø 30 mm, H 115 mm, Max. 17000xg	300 pezzi



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

