

**Dosatore per bottiglia Varispenser® 2x****Eppendorf AG**

Dispensatore per bottiglie di reagenti. Dispensazione sicura e semplice di liquidi da flaconi. Varispenser 2 e Varispenser 2x sono la scelta ideale per la dispensazione di aliquote di liquido da flaconi di alimentazione e di reagente. La tecnologia avanzata del bordo di tenuta permette la dispensazione di quasi tutti i tipi di acidi, soluzioni o basi. Disponibili in 6 misure per dispensare i volumi tra 0,2 # 100 mL, questi apparecchi sono progettati per dispensare costantemente i volumi selezionati senza perdita del reagente. Tutti i dosatori Varispenser sono completamente autoclavabili, per una massima sicurezza nel vostro laboratorio. Varispenser 2x ha una valvola di dispensazione di ritorno che evita la perdita del reagente durante la ventilazione.



Code	Description	Packaging
LLG06272404	Dispensatore Varispense 2 plus 1-10ml Adatt.GL25, GL28/S 28, GL32, GL38, S40	1 pz
LLG06272409	Dispensatore Varispencer 2 1-10ml Adatt.GL25, GL28/S 28, GL32, GL38, S40	1 pz
LLG06272405	Dispensatore Varispenser® 2x, 2,5-25 ml, Adattatore GL 32, GL 38, S 40	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

