

**COD-Cuvette test standard, 16mm****AQUALYTIC**

COD-Cuvette test standard, 16mm. Adatto per vari fotometri, grazie alle proprietà chimiche (es. Aqualytic, Hach) con provette diametro 16mm.



| Code        | Description                                                                                                    | Packaging |
|-------------|----------------------------------------------------------------------------------------------------------------|-----------|
| LLG09950533 | COD-Cuvette, 0-150 mg/L O <sub>2</sub>                                                                         | 25 pz     |
| LLG09950534 | COD-Cuvette, 0-1.500 mg/L                                                                                      | 25 pz     |
| LLG09950535 | COD-Cuvette, 0-15.000 mg/L O <sub>2</sub>                                                                      | 25 pz     |
| LLG06292149 | Cuvetta per test COD, 16 mm, Tipo CSB LMR, Range di misura 15 ... 300 mg/l O <sub>2</sub> , UN 2922, 8, II (E) | 25 tests  |
| LLG06242606 | Test in cuvetta CSB/COD Vario, 0-1500 mg/l O <sub>2</sub> , 150 test, UN 2922, 8, II, (E), 16 mm               | 150 pz    |
| LLG06242848 | Cuvette COD test standard, 16 mm, Tipo CSB MR senza mercurio, Range di misura 0 ... 1500 mg/l O <sub>2</sub>   | 25 pz     |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

