

**Elettrodi pH InLab® Expert Pro****Mettler-Toledo Online**

Elettrodi pH InLab® Expert Pro. Elettrodi combinati pH con stelo in PEEK chimicamente resistente, sistema di riferimento ARGENTAL#. Cavo fisso BNC/RCA (Cinch). InLab® Expert Pro ISM adatto per SevenCompact e SevenExcellence.



Code	Description	Packaging
LLG09041504	Elettrodo inlab expert ntc30	1 pz.
LLG09041506	Elettrodo pH InLab Expert Pro	1 pz.
LLG09041508	Elettrodo pH InLab Expert Pro ISM	1 pz.



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

