



## Telecamera per microscopio Moticom 1080

Motic

Telecamera per microscopio Moticom 1080, risoluzione 2,0 MP, formato sensore 1/2,8".



| Code        | Description                                                            | Packaging |
|-------------|------------------------------------------------------------------------|-----------|
| LLG06282488 | Telecamera per microscopio Moticam 1080, risol. 2,0 MP, sensore 1/2,8" | 1 pz      |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

