



Albero agitatore KPG

Albero agitatore KPG. DURAN®. Intercambiabile.



Code	Description	Packaging
LLG09197311	Asta d'agitazione KPG, DURAN, Tipo WS 10, Lungh. 320 mm, per Ø collo 25 mm	1 pz
LLG09197312	Asta d'agitazione KPG, DURAN, Tipo WS 10, Lungh. 370 mm, per Ø collo 25 mm	1 pz
LLG09197313	Asta d'agitazione KPG, DURAN, Tipo WS 10, Lungh. 410 mm, per Ø collo 25 mm	1 pz
LLG09197314	Asta d'agitazione KPG, DURAN, Tipo WS 10, Lungh. 440 mm, per Ø collo 25 mm	1 pz
LLG09197321	Asta d'agitazione KPG, DURAN, Tipo WG 10, Lungh. 320 mm, per Ø collo 60 mm	1 pz
LLG09197322	Asta d'agitazione KPG, DURAN, Tipo WG 10, Lungh. 370 mm, per Ø collo 60 mm	1 pz



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

