

**Palloni Kjeldahl****Isolab**

Palloni Kjeldahl. Prodotto in vetro borosilicato 3.3. I colli NS sono in conformità con gli standard DIN 12242.



| Code        | Description                              | Packaging |
|-------------|------------------------------------------|-----------|
| LLG09140385 | Palloni Kjeldahl, 100 ml, cono 19/26 NS  | 1 pz.     |
| LLG09140389 | Palloni Kjeldahl, 1000 ml, cono 29/32 NS | 1 pz.     |
| LLG09140386 | Palloni Kjeldahl, 250 ml, cono 29/32 NS  | 1 pz.     |
| LLG09140387 | Palloni Kjeldahl, 500 ml, cono 29/32 NS  | 1 pz.     |
| LLG09140388 | Palloni Kjeldahl, 750 ml, cono 29/32 NS  | 1 pz.     |



**CARLO ERBA Reagents S.r.l.**  
 VIA R. MERENDI, 22  
 20007 CORNAREDO (MI)  
 TEL. +39 02 93 99 190  
 FAX +39 02 93 991 001  
[www.carloerbareagents.com](http://www.carloerbareagents.com)



**CARLO ERBA**  
*Reagents operates with  
 a Certified Quality  
 Management System*

