

**Test Visocolor Alpha****Macherey Nagel**

I test VISOCOLOR® alpha possono essere di tipo colorimetrico o titrimetrico. Facili e sicuri da maneggiare, in quanto, spesso, è necessario un solo reagente per l'esecuzione del test. Le bottiglie di reagente sono confezionate in pratici blister. Il retro del cartone slitta per chiusura/apertura della confezione e fornisce tutte le informazioni necessarie per il test: istruzioni in 6 lingue con pittogrammi, scala di comparazione per valutazione. La confezione blister è dotata di foro per facilitare lo stoccaggio.



Code	Description	Packaging
H2O935012	VISOCOLOR ALPHA Ammonio - 0/0.2/0.5/1/2/3 mg/L NH ₄ ⁺	1 kit
H2O935016	VISOCOLOR ALPHA Alcalinità carbonato - 1 goccia = 1°d (corrisponde a 18mg/L CaCO ₃)	1 kit
H2O935019	VISOCOLOR ALPHA Cloro - 0.25/0.5/1.0/1.5/2.0 mg/L Cl ₂	1 kit
H2O935042	VISOCOLOR ALPHA Durezza totale - 1 goccia=1.3°e (corrisponde a 18mg/L CaCO ₃)	1 kit
H2O935065	VISOCOLOR ALPHA Nitrato - 2/8/15/30/50 mg/L NO ₃ ⁻	1 kit
H2O935066	VISOCOLOR ALPHA Nitriti - 0.05/0.10/0.25/0.5/1.0 mg/L NO ₂ ⁻	1 kit
H2O935075	VISOCOLOR ALPHA pH 5-9	1 kit
H2O935079	VISOCOLOR ALPHA Fosfato - 2/5/10/15/20 mg/L PO ₄ ³⁻	1 kit
H2O935080	VISOCOLOR ALPHA Durezza residua - 0.00/0.04/0.08/0.15/0.30 °d	1 kit



CARLO ERBA Reagents S.r.l.
 VIA R. MERENDI, 22
 20007 CORNAREDO (MI)
 TEL. +39 02 93 99 190
 FAX +39 02 93 991 001
www.carloerbareagents.com



CARLO ERBA
*Reagents operates with
 a Certified Quality
 Management System*

